



Behavior Hub

Lab quick start

1.2.6-r2 - LAB QUICK START

Before an experiment: make a short test recording and check the saved video for clear focus, framing and lighting.

1. Open the app and a project

Extract the complete ZIP. Open **BehaviorHub.exe**, then choose **New Project** or **Open Project**. Keep the EXE and _internal together.

2. Set duration and video names

Choose behavior and duration. For **Custom**, enter the behavioral test name. MP4 names use **date + time + behavioral test + Mouse ID**. Blank Mouse ID receives **Auto_ID01**, **Auto_ID02**, **Auto_ID03**. Auto_ID means an automatically assigned mouse ID when the field is blank. Numbering is saved with the project and continues after restart. Test mode uses seconds.

3. Check camera and focus

Select the camera, resolution and target FPS. Click **Open preview**. Choose **Manual** and drag the focus slider until clear. Changes apply and save automatically. **Zoom and flip** offers 1-4x zoom, pan directions and two flip options. Framing updates instantly without closing preview. **Reset** restores the full view. Both preview and saved MKV/MP4 use the adjusted view; no uncropped copy is kept. Controls lock during recording.

4. Record

Click **Start recording**. The timer stops recording automatically. Use **Stop early and save** only when needed. Names are fixed at recording start; duplicates get a number.

5. Process and review

Close preview or click **Process pending videos**. Wait for **MP4: verified**. Check any **REVIEW** result before analysis.

6. Find all new MP4 files

Open the project **MP4** folder. Example: **20261002_143025_OFT_Mouse01.mp4**. TEST videos have a **TEST_** prefix. Originals and timing reports stay in the session folders. Older videos keep their original locations.

7. Analyze and back up

For DLC, select one behavior and exclude TEST videos. Save matching DLC CSVs beside MP4s. Copy the **whole project folder** for backup. Use **Exit application** when finished.